EN Sterilization wraps

4501 - Sterilization wrap made of nonwoven SMS polypropylene, blue

4502 - Sterilization wrap made of nonwoven SMS polypropylene, green (alcohol repellent)

4503 - Sterilization wraps made of nonwoven SMS polypropylene, blue and green interleaved

450 - Creped medical grade paper (white, blue or green)

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450 - Papel crepado de grau médico (branco, azul ou verde)

EN Sterilization wraps

INTENDED PURPOSE AND CLINICAL BENEFIT

Sterilization wraps are used for packaging of medical devices for sterilization in order to maintain the sterility of terminally sterilized medical devices until the point of use, in order to prevent patient infection.

Nonwoven SMS polypropylene sterilization wraps are designed for sterilization by steam (suitable for 121°C or 134°C; process in conformance with ISO 17665), ethylene oxide (process in conformance with ISO 11135) and vaporized hydrogen peroxide - VH2O2 (process in conformance with ISO 22441).

Crepe paper is designed for sterilization by steam (suitable for 121°C or 134°C; process in conformance with ISO 17665), ethylene oxide (process in conformance with ISO 11135) and radiation (process in conformance with ISO 11137).

PRODUCT INFORMATION

In compliance with the requirements from ISO 11607 and EN 868-2, aiming to ensure microbial barrier properties and the subsequent sterility maintenance. Designed for an easy and versatile aseptic presentation of the sterilized medical devices.

INSTRUCTIONS FOR USE

1- Visually inspect the sterilization wraps before use, in order to confirm there are no defects or foreign elements.
2- Select the size of the sterilization wrap according to size, shape and weight of the medical devices to be wrapped or based on guidelines within the organization.

The size should be selected to achieve adequate coverage of the articles being wrapped but it should not be so big that it has to be wrapped several times around the medical devices, as this could difficult the penetration of the sterilant agent.

3- Perform the proper wrapping technique, in order to allow aseptic presentation of the sterile content. The most common examples of wrapping techniques are the envelope fold and the square fold.

When two layers of sterilization wraps are used, the simultaneous double-wrapping technique (that refers to wrapping two layers of wrap material together simultaneously), or the sequential wrapping technique (that refers to wrapping two layers of wrap material individually) can be performed. Note that it is essential to wrap the medical devices securely to prevent gaps, billowing and air pockets, but do not wrap too tightly as this could create holes or tears in the wrap. Furthermore, the sterilization wrap should be large enough to accommodate movement of the wrap during the sterilization cycle without ripping or tearing.

4- Close the wrapped packages with an appropriate closure system. The most common closure system is the indicator tape.

Note that the closure system should maintain the integrity of the package.

Do not use closures systems that avoid the penetration of the sterilant agent nor that compresses the package (e.g., ropes, strings, elastic bands, paperclips, staples or similar items).

5- Submit the package to the appropriate validated sterilization method.

6- After sterilization, visually inspect the package to ensure it is intact and the product is clean and dry. Wet packages are considered non-sterile.

7- Handle and storage the product carefully in order to maintain the integrity of the packaging.

If sterilization is carried out by a subcontracted unit, it is recommended that the packages are protected from possible contamination with a sterility maintenance cover.

WARNINGS AND PRECAUTIONS

Sterilization wraps are intended to be used exclusively by sterilization professionals for final presentation of devices to healthcare professionals.

Articles to be packed must be completely clean and dry. Avoid sharp edges and other surfaces that might disrupt the package. If it is not possible to avoid them, they should be protected.

Keep the sterilization wraps in a clean and dry environment. We decline any responsibility in case of:

- inappropriate use;
- inappropriate storage or handling conditions;
- inappropriate sterilization conditions or sterilization method.

Do not re-use. Sterilization wraps are intended to be used one single time to prevent risks of compromised microbial barrier properties and the subsequent compromised sterility.

Do not re-sterilize.

Visually inspect the sterile barrier system and the device prior to use. Do not use the device if you detect any breach in packaging integrity or any sign of the device deterioration that may compromise its use and/or functionality.

Sterilization wraps maintain the product sterility up to 180 days.

Dispose according to local legal regulations for urban waste.

The occurrence of serious incidents is not expected with the use of sterilization wraps, but we must inform that serious incidents related with the use of the device, that are life-threatening or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your country.

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