
Infusion Adapter

Infusion Adapter (C100-0)

| [REF 515078](#)

Infusion Adapter (C100-0 MULTI)

| [REF 515079](#)

[en] English

DEVICE DESCRIPTION: The BD PhaSeal™ Optima Infusion Adapter is sterile, single-use, and part of the BD PhaSeal™ Optima System.

INTENDED PURPOSE: The Infusion Adapter is an IV container spike that attaches to an IV line and is intended for the closed preparation and/or administration of hazardous parenteral drugs when used with a BD PhaSeal™ Optima Injector.

INTENDED USER: The Infusion Adapter is intended for use by healthcare workers who prepare and administer hazardous parenteral drugs.

INTENDED PATIENT POPULATION: Adults and pediatrics.

PATIENT SELECTION CRITERIA: Patients who receive parenteral hazardous drugs.

INTENDED USE / INDICATIONS FOR USE:

The BD PhaSeal™ Optima System is an airtight and leakproof closed system drug transfer device (CSTD) that mechanically prohibits the transfer of environmental contaminants into the system and the escape of drug vapor concentrations outside the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols and spills. The BD PhaSeal™ Optima System also prevents microbial ingress for up to 168 hours.

PRECAUTIONS & WARNINGS:

- Do not use if package is damaged.
- Do not use if protection caps are displaced.
- Do not reuse to avoid contamination.
- Follow standards of practice and facility policy on sterile and hazardous drug compounding which includes but is not limited to the use of a primary engineering control that provides an ISO Class 5 environment or better.
- The Infusion Adapter is compatible with standard infusion bags/containers and the spike is designed according to relevant ISO standards.

CONTRAINDICATIONS: Do not use with drugs contraindicated for polypropylene (PP). Not for fat emulsions or blood components.

GENERAL GUIDELINES:

A non-vented IV set should be used together with the Infusion Adapter. If a vented IV set is used, make sure that the air inlet is closed throughout the infusion procedure. The Infusion Adapter may be used for up to 10 penetrations.

EU ONLY: Users should report any serious incident related to the device to the Manufacturer and National Competent Authority.

Note: The device is provided sterile. Follow your institutional guidelines for disinfecting prior to initial entry or reconnection. Disinfect the device by scrubbing with 70% isopropyl alcohol pad and by wiping with friction and using a twisting motion for at least three seconds. Allow wiped surfaces to dry for a minimum of one minute.

Component Materials - PP, TPE, Isoprene

INSTRUCTIONS FOR USE:

- Remove the protection cap from the spike. Push the spike fully into the IV container.
- Remove the cap from the Infusion Adapter connection interface. Slide the injector straight onto the connection interface until the two pieces snap together.
- Perform drug transfer.
- Pull the injector straight back until the two pieces are disconnected.
- Open the blue protection lid on the Infusion Adapter. While holding the bag spike port, connect the IV set to the Infusion Adapter with the IV container in horizontal position.
- Dispose without removing the components in appropriate container after single use.

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