

Parallel jaw closure. **Even tissue compression.**

PROXIMATE® Reloadable Linear Stapler (TX)



Consistent staple formation

- Parallel jaw closure helps to compress tissue evenly within the jaws of the instrument for consistent staple formation

Easy to use

- Intermediate locking position facilitates proper placement of stapler on tissue
- One handed operation to close, fire and release

Applicable specialties

- General
- Colorectal
- Thoracic

Product Overview

- Two staggered rows of titanium alloy staples in 30mm and 60mm sizes.
- Three different reloads accommodates a variety of tissue thicknesses

Device Specifications

Description	Device Code	Tissue	Staple rows	Staple line length	Number of staples	Max firings	Qty/box
Proximate Linear Stapler 30mm Reload-Vascular	TX30V	<input type="checkbox"/> Vascular	3	30mm	23	8	3
Proximate Linear Stapler 30mm	TX30B	<input checked="" type="checkbox"/> Regular	2	30mm	11	8	3
	TX30G	<input checked="" type="checkbox"/> Thick					
Proximate Linear Stapler 60mm	TX60B	<input checked="" type="checkbox"/> Regular	2	60mm	21	8	3
	TX60G	<input checked="" type="checkbox"/> Thick					



Reload Specifications

Description	Reload code	Tissue	Closed staple height (mm)	Wire diameter	Open staple height	Qty/box
Proximate Linear Stapler 30mm Reload-Vascular	XR30V	<input type="checkbox"/> Vascular	1.0mm	0.20mm	2.4mm	12
Proximate Linear Stapler 30mm Reload	XR30B	<input checked="" type="checkbox"/> Regular	1.5mm	0.23mm	3.6mm	12
	XR30G	<input checked="" type="checkbox"/> Thick	2.0mm	0.27mm	4.7mm	
Proximate Linear Stapler 60mm Reload	XR60B	<input checked="" type="checkbox"/> Regular	1.5mm	0.23mm	3.6mm	12
	XR60G	<input checked="" type="checkbox"/> Thick	2.0mm	0.27mm	4.7mm	



How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway: www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX): www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI): JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer Support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at customersupport@eesus.jnj.com or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.