

# VetO2+ Silent

## Oxygen Concentrator



**EVEN MORE  
EFFECTIVE,  
QUIETER,  
SMALLER!**

# VetO2+ Silent

## Oxygen concentrator



### When to use an oxygen concentrator?

In your clinic, the oxygen concentrator is part of the important equipment. Much more advantageous than oxygen cylinders, it can be used:

- ▶ Either to bring the oxygen that the animal needs, via a mask or a tracheal intubation or in an oxygen cage for oxygen therapy.
- ▶ Or during an anesthesia, supplying oxygen as vector gas

In both situations, the oxygen concentrator has considerable advantages over the use of oxygen cylinders:

- √ No risk of explosion
- √ No need to stock or purchase often oxygen cylinders

The oxygen concentrator produces oxygen concentrated up to 90-95% anytime. It is a reliable, robust and easy-to-use equipment.



### Why choose VetO2+ Silent oxygen concentrator?

- ▶ Quiet: only 40 dB!
- ▶ Equipped with a humidification flask including a pressure relief valve
- ▶ Compatible with most anesthesia systems
- ▶ Oxygen concentration management, electronic hours meter, timer programming up to 40 hours



### VetO2+ Silent technical features:

Flow:	Up to 5l/min
Oxygen concentration:	93% ±3
Maximum outlet pressure:	0.04-0.07 Mpa
Power supply:	AC220V/110V ±10%
Sound level:	40 dB max.
Weight:	16 Kg
Dimensions:	35× 28 × 51 cm
Classification:	II Class Type B
LCD display:	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10 min. to 40 hours
Alarm:	Power failure, high and low pressure, low purity



15, rue du Noroît  
22100 Taden - France  
Tel. +33 (0)2 96 85 86 76

[welcome@manomedical.com](mailto:welcome@manomedical.com)

[www.manomedical.com](http://www.manomedical.com)