MINIDOP ES-100VX

POCKET DOPPLER
WITH GREAT
SOUNDS





For Vascular and Fetal

- HAND-HELD DESIGN
- GREAT SENSITIVITY
- WIDE PROBE SELECTION
 For user flexibility. (2, 4, 5, 8 and 10 MHz)
- CONVENIENT PROBE ACTIVATION BUTTON
 Turns ES-100VX ON and OFF.
- AUTOMATIC POWER "OFF"

PROBE SELECTION



CLINICAL APPLICATIONS

- » Peripheral vascular procedures
- » Venous compressions
- » Fetal heart rate

- » Blood Pressure segmental studies
- » Penile & Digit systolic pressures
- » Flow detection in recovery room

SPECIFICATIONS

Power: DC9V, A 9 volt alkaline

battery

Battery life: Approx. 3 hours

Automatic shut-off

Multi frequency: 2, 4, 5, 8, 10 MHz
Probe button: Power ON/OFF
Speaker output: 200 mW or more

Outputs (3.5mm jack)

Headset: Speaker cut off.

Dimensions: 75(W)x140(D)x25(H) mm

Weight: Approx. 270 grams

(including probe)

Accessories: Carrying case, camera

(optional) strap, gel, headset, battery



