

veterinary technology for life **Eickeneyer**®

For Your Eyes Only



AlgerBrush II Corneal Burr

The AlgerBrush II is a battery powered, handheld and low torque motorised burr. It has been designed to remove foreign bodies from the cornea or sclera of patients.

The burr in this rotary instrument is used as a "brush" rather than a drill and the foreign bodies are thus "brushed" from the eye leaving a smooth surface, which in turn heals much faster than with other instruments that could be used for this purpose.

According to the paper "Veterinary Ophthalmology" [(2013) 16, 2, 83-88], it has been successfully used to treat spontaneous corneal epithelial defects (SCCEDs) in dogs in combination with a bandage contact lens.

Includes 3.5 mm round, medium grit diamond burr

170010

AlgerBrush II Chuck and Burr

3.5 mm round, medium grit diamond burr

170011

AlgerBrush II Corneal Burrs

1 mm

5 per pack

170012

AlgerBrush II Chuck and Burr2.5 mm round, medium grit diamond burr

170013

AlgerBrush II Diamond Grinding Wheel (4 mm) Fine

Convenient for recurrent corneal erosion

170014

EICKEMEYER[®] Eltastraße 8 78532 Tuttlingen Germany

T +49 7461 96580 0 F +49 7461 96580 90 export@eickemeyer.com