

Gauze balls



- ✓ Cotton, with elastic ring.
- ✓ For wound dressing procedures, skin disinfection procedures, swabbing and absorbing during surgery.
- ✓ X-Ray detectable thread.
- ✓ Raw edges of the fabric are turned inside and kept with the help of an elastic ring made of silicone.
- ✓ 17 threads/sq.cm with approx. 23g/sq.m
20 threads/sq.cm with approx. 27g/sq.m



REF	Th./sq.cm	REF	Weight g/pc.	Diameter (mm)		X-Ray	Gauze size	Pcs./Bag	Pcs./Tr.carton
				X	Y				
429-304	17	D	0,37	25	35		180cm ²	1.000	2.000
429-305	17	E	0,71	30	40		340cm ²	750	1.500
429-306	17	F	0,83	40	50		400cm ²	750	1.500
429-307	17	G	1,63	50	60		780cm ²	500	1.500
429-308	17	H	1,87	55	70		900cm ²	400	1.200
429-309	17	I	2,92	65	80		1400cm ²	250	1.000
429-310	17	J	3,32	75	90		1600cm ²	200	800
429-311	17	K	3,78	85	100		1800cm ²	200	800
429-315	17	D	0,37	25	35	X	180cm ²	1.000	2.000
429-316	17	E	0,71	30	40	X	340cm ²	750	1.500
429-317	17	F	0,83	40	50	X	400cm ²	750	1.500
429-318	17	G	1,63	50	60	X	780cm ²	500	1.500
429-319	17	H	1,87	55	70	X	900cm ²	400	1.200
429-326	20	D	0,44	25	35		180cm ²	1.000	2.000
429-327	20	E	0,83	30	40		340cm ²	750	1.500
429-328	20	F	0,98	40	50		400cm ²	750	1.500
429-329	20	G	1,91	50	60		780cm ²	500	1.500
429-330	20	H	2,2	55	70		900cm ²	400	1.200
429-331	20	I	3,42	65	80		1400cm ²	250	1.000
429-332	20	J	3,9	75	90		1600cm ²	200	800
429-333	20	K	4,44	85	100		1800cm ²	200	800
429-337	20	D	0,44	25	35	X	180cm ²	1.000	2.000
429-338	20	E	0,83	30	40	X	340cm ²	750	1.500
429-339	20	F	0,98	40	50	X	400cm ²	750	1.500
429-340	20	G	1,91	50	60	X	780cm ²	500	1.500
429-341	20	H	2,2	55	70	X	900cm ²	400	1.200
429-342	20	I	3,42	65	80	X	1400cm ²	250	1.000
429-343	20	J	3,9	75	90	X	1600cm ²	200	800
429-344	20	K	4,44	85	100	X	1800cm ²	200	800